

ECT Step by Step

Lab Scales

7/24/08

Lab Scales

For New Test

Select FGIS-904 Lab Scale Test from Dropdown next to Create New Test

To finish an Open test use Search or enter the Test Reference No.

USDA United States Department of Agriculture
Grain Inspection, Packers & Stockyards Administration

Staging

Home | About GIPSA | Help | Contact Us | Logoff

You are here: ECT Menu 1.0.5

Checktest Home

Create New Checktest

Search For Existing Checktest

Search By Test Reference No.

Equipment

Reports

E-mail Distribution Lists

Administration

Search

Search for Checktest

Browse by Subject

- Equipment
- Reports
- E-mail Distribution Lists
- Administration

FGIS-904 Laboratory Scale Test

FGIS-923 Moisture Meter Test

FGIS-924 Barley Pearler

FGIS-924 Sieve Test

FGIS-925 Rice Checktest Form

FGIS-927 Test Weight Check Test

FGIS-928 Dockage Check Test

FGIS-936 Sampler Condition Report

FGIS-965 Scale Test Report (Hopper)

FGIS-965-1 Scale Test Report (Railroad Track)

FGIS-965-2 Scale Test Report (Vehicle)

FGIS-Falling Number

FGIS-Mycotoxin

FGIS-NIRT Daily Barley

FGIS-NIRT Daily Corn

FGIS-NIRT Daily Soybean

FGIS-NIRT Daily Wheat

GIPSA | USDA.gov
FOIA | Accessibility Statement | Privacy Policy | Non-Discrimination Statement | Information Quality | FirstGov | White House

15 in the mail. RTDRA | USDA.gov

FGIS-904 Laboratory Scale Test General Information

Test Reference Number:

Test Date:

Service Point

FO OA LOCATION

Laboratory Scale

Capacity g X Division Size 1 (d) g

Dual Range Scales Only:

Break Point g Division Size 2 (d) g

Class II

Class II	
No. of Div.	Tol.
0-5,000	1d
5,001-20,000	2d
20,001+	3d

 Class III

Class III	
No. of Div.	Tol.
0-500	1d
501-2,000	2d
2,001-4,000	3d
4,001+	5d

Type	d	Test Loads Required	Tol.
<input type="radio"/> Precision, etc.	0.01 g	1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 20, 50, 100 g	0.02 g
<input type="radio"/> Moisture, etc.	0.1	above plus 150, 200, 250, 300, 400 g	0.2g
<input type="radio"/> General, etc.	1	above plus 500, 600, 700, 1000, 1500 g	1

Blank General Page

Steps to Complete:

1. Enter Test Date
2. Enter Service Point
3. Click on Add Equipment
4. Use Equipment search to find equipment and then Select Equipment
5. Select Capacity from Dropdown
6. Select Division Size from Dropdown
7. Only Dual range scales require use of breakpoint and Division Size 2
8. Select either Class II, Class III, Precision, Moisture, General
Only select one.
9. Click on Save
10. Review page
11. Click on Next

FGIS-904 Laboratory Scale Test General Information

Test Reference Number: 14043

Checktest record saved

Test Date:

Service Point

FO **OA** **LOCATION**
 FGIS - Wichita Field Office KANSAS CITY, MO

Laboratory Scale

Select	SP Code	Serial #	Make	Model
[Remove]	244121	G-31359	Mettler	PM400/49

Capacity g X Division Size 1 (d) g

Dual Range Scales Only:

Break Point g Division Size 2 (d) g

<input checked="" type="radio"/> Class II	<table border="1"> <thead> <tr> <th colspan="2">Class II</th> </tr> <tr> <th>No. of Div.</th> <th>Tol.</th> </tr> </thead> <tbody> <tr> <td>0-5,000</td> <td>1d</td> </tr> <tr> <td>5,001-20,000</td> <td>2d</td> </tr> <tr> <td>20,001+</td> <td>3d</td> </tr> </tbody> </table>	Class II		No. of Div.	Tol.	0-5,000	1d	5,001-20,000	2d	20,001+	3d	<input type="radio"/> Class III	<table border="1"> <thead> <tr> <th colspan="2">Class III</th> </tr> <tr> <th>No. of Div.</th> <th>Tol.</th> </tr> </thead> <tbody> <tr> <td>0-500</td> <td>1d</td> </tr> <tr> <td>501-2,000</td> <td>2d</td> </tr> <tr> <td>2,001-4,000</td> <td>3d</td> </tr> <tr> <td>4,001+</td> <td>5d</td> </tr> </tbody> </table>	Class III		No. of Div.	Tol.	0-500	1d	501-2,000	2d	2,001-4,000	3d	4,001+	5d
Class II																									
No. of Div.	Tol.																								
0-5,000	1d																								
5,001-20,000	2d																								
20,001+	3d																								
Class III																									
No. of Div.	Tol.																								
0-500	1d																								
501-2,000	2d																								
2,001-4,000	3d																								
4,001+	5d																								

Type	d	Test Loads Required	Tol.
<input type="radio"/> Precision, etc.	0.01 g	1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 20, 50, 100 g	0.02 g
<input type="radio"/> Moisture, etc.	0.1	above plus 150, 200, 250, 300, 400 g	0.2g
<input type="radio"/> General, etc.	1	above plus 500, 600, 700, 1000, 1500 g	1

Last change made on 7/2/2008 11:45:10 AM by user Mary Vick.

Completed General Page


FGIS-904 Laboratory Scale Test Test Information

General Test

Test Reference Number: 14043

Blank Test Page

1. Enter date test weights last checked.
2. Record Yes (No will fail test) for zero balancing the scale. Click over correct radio button.
3. Unless this is a mechanical scale select N/A by clicking over radio button. If mechanical Yes is required to pass.
4. Enter Shift test load
5. Enter Actual Test values.
6. Precision of results must match the division size
7. Click on Save
8. Review
9. Click on Submit

1. **Class F weights** used in this exam were last tested on this date 
2. **Zero balance** the scale before the sensitivity test, shift test and increasing load test.
[Yes No]
3. **Sensitivity** Mechanical scales only. At zero, add weight equal to one scale division. The scale indicator should move to the top of the trig loop (or readings should change by at least one division). Repeat the test at the scale's maximum test load.
[Yes No N/A]

Shift Test				
Shift Load Test: <input type="text"/> (used as expected shift test reading)				
	Quadrant 1		Quadrant 2	
Reading	<input type="text"/>		<input type="text"/>	<input type="text"/>

Increasing Load							
Test Load	Indication	Test Load	Indication	Test Load	Indication	Test Load	Indication
1g	<input type="text"/>	7g	<input type="text"/>	100g	<input type="text"/>	500	<input type="text"/>
2	<input type="text"/>	8	<input type="text"/>	150	<input type="text"/>	600	<input type="text"/>
3	<input type="text"/>	9	<input type="text"/>	200	<input type="text"/>	700	<input type="text"/>
4	<input type="text"/>	10	<input type="text"/>	250	<input type="text"/>	1000	<input type="text"/>
5	<input type="text"/>	20	<input type="text"/>	300	<input type="text"/>	1500	<input type="text"/>
6	<input type="text"/>	50	<input type="text"/>	400	<input type="text"/>		<input type="text"/>

Decreasing Load	
Test Load	Indication
400	<input type="text"/>
200	<input type="text"/>
50	<input type="text"/>

Balance Change	
Test Load	Indication
0	<input type="text"/>


Approved 	Disapproved 	Tested By: <input type="text"/>	
--	---	---------------------------------	--

FGIS-904 Laboratory Scale Test Test Information

General | **Test**

Test Reference Number: 14043

Indication for increasing test load 5 does not show the precision of exactly 2 digit(s) to the right of the decimal point
Indication for balance test load 0 does not show the precision of exactly 2 digit(s) to the right of the decimal point
Results by (inspector) license code required

1. **Class F weights** used in this exam were last tested on this date 

2. **Zero balance** the scale before the sensitivity test, shift test and increasing load test.
[Yes No]

3. **Sensitivity** Mechanical scales only. At zero, add weight equal to one scale division. The scale indicator should move to the top of the trig loop (or readings should change by at least one division). Repeat the test at the scale's maximum test load.
[Yes No N/A]

Shift Test				
Shift Load Test: <input type="text" value="200.00"/> (used as expected shift test reading)				
	Quadrant 1	Quadrant 2	Quadrant 3	Quadrant 4
Reading	<input type="text" value="200.00"/>	<input type="text" value="200.01"/>	<input type="text" value="200.00"/>	<input type="text" value="200.00"/>

Increasing Load							
Test Load	Indication	Test Load	Indication	Test Load	Indication	Test Load	Indication
1g	<input type="text" value="1.00"/>	7g	<input type="text" value="7.00"/>	100g	<input type="text" value="100.00"/>	500	<input type="text"/>
2	<input type="text" value="2.00"/>	8	<input type="text" value="8.00"/>	150	<input type="text" value="150.00"/>	600	<input type="text"/>
3	<input type="text" value="3.00"/>	9	<input type="text" value="9.00"/>	200	<input type="text" value="200.01"/>	700	<input type="text"/>
4	<input type="text" value="4.00"/>	10	<input type="text" value="10.00"/>	250	<input type="text" value="250.00"/>	1000	<input type="text"/>
5	<input type="text" value="5.0"/>	20	<input type="text" value="20.00"/>	300	<input type="text" value="300.00"/>	1500	<input type="text"/>
6	<input type="text" value="6.00"/>	50	<input type="text" value="50.00"/>	400	<input type="text" value="399.99"/>		<input type="text"/>

Error Messages

The system will check for missing information and will check the data for correct precision.

It will not let you submit results until these errors are corrected.

Test that passed

Click on Cancel to go back to Home page or use breadcrumb menu at the top

Note: Approved radio button is now selected.

You are here: ECT Menu > FGIS-904 Laboratory Scale Form 1.0.5

FGIS-904 Laboratory Scale Test Test Information

Test Reference Number: 14043

Checktest submitted and APPROVED

- Class F weights** used in this exam were last tested on this date
- Zero balance** the scale before the sensitivity test, shift test and increasing load test.
 Yes No]
- Sensitivity** Mechanical scales only. At zero, add weight equal to one scale division. The scale indicator should move to the top of the trig loop (or readings should change by at least one division). Repeat the test at the scale's maximum test load.
 Yes No N/A]

Shift Test				
Shift Load Test: <input type="text" value="200.00"/> (used as expected shift test reading)				
	Quadrant 1	Quadrant 2	Quadrant 3	Quadrant 4
Reading	<input type="text" value="200.00"/>	<input type="text" value="200.01"/>	<input type="text" value="200.00"/>	<input type="text" value="200.00"/>

Increasing Load							
Test Load	Indication	Test Load	Indication	Test Load	Indication	Test Load	Indication
1g	<input type="text" value="1.00"/>	7g	<input type="text" value="7.00"/>	100g	<input type="text" value="100.00"/>	500	<input type="text"/>
2	<input type="text" value="2.00"/>	8	<input type="text" value="8.00"/>	150	<input type="text" value="150.00"/>	600	<input type="text"/>
3	<input type="text" value="3.00"/>	9	<input type="text" value="9.00"/>	200	<input type="text" value="200.01"/>	700	<input type="text"/>
4	<input type="text" value="4.00"/>	10	<input type="text" value="10.00"/>	250	<input type="text" value="250.00"/>	1000	<input type="text"/>
5	<input type="text" value="5.00"/>	20	<input type="text" value="20.00"/>	300	<input type="text" value="300.00"/>	1500	<input type="text"/>
6	<input type="text" value="6.00"/>	50	<input type="text" value="50.00"/>	400	<input type="text" value="399.99"/>		<input type="text"/>

Decreasing Load	
Test Load	Indication
400	<input type="text" value="399.99"/>
200	<input type="text" value="200.00"/>
50	<input type="text" value="50.00"/>

Balance Change	
Test Load	Indication
0	<input type="text" value="0.00"/>

Approved
 Disapproved
 Tested By:

FGIS-904 Laboratory Scale Test Test Information

General

Test

Test Reference Number: 14022

Indication for increasing test load 100 outside of tolerance range of +/- 0.020
Checktest submitted and DISAPPROVED

1. **Class F weights** used in this exam were last tested on this date

2. **Zero balance** the scale before the sensitivity test, shift test and increasing load test.
[Yes No]

3. **Sensitivity** Mechanical scales only. At zero, add weight equal to one scale division. The scale indicator should move to the top of the trig loop (or readings should change by at least one division). Repeat the test at the scale's maximum test load.
[Yes No N/A]

Test that Failed

Information show in red is not saved as part of the record.

I usually copy that and paste into remarks as a reminder of why a test failed.

Shift Test				
Shift Load Test: <input type="text" value="55.00"/> (used as expected shift test reading)				
	Quadrant 1		Quadrant 2	
Reading	<input type="text" value="55.00"/>		<input type="text" value="55.00"/>	<input type="text" value="55.00"/>

Increasing Load							
Test Load	Indication	Test Load	Indication	Test Load	Indication	Test Load	Indication
1g	<input type="text" value="1.00"/>	7g	<input type="text" value="7.00"/>	100g	<input type="text" value="100.05"/>	500	<input type="text"/>
2	<input type="text" value="2.00"/>	8	<input type="text" value="8.00"/>	150	<input type="text"/>	600	<input type="text"/>
3	<input type="text" value="3.00"/>	9	<input type="text" value="9.00"/>	200	<input type="text"/>	700	<input type="text"/>
4	<input type="text" value="4.00"/>	10	<input type="text" value="10.00"/>	250	<input type="text"/>	1000	<input type="text"/>
5	<input type="text" value="5.00"/>	20	<input type="text" value="20.00"/>	300	<input type="text"/>	1500	<input type="text"/>
6	<input type="text" value="6.00"/>	50	<input type="text" value="50.02"/>	400	<input type="text"/>		<input type="text"/>

Decreasing Load	
Test Load	Indication
400	<input type="text"/>
200	<input type="text"/>
50	<input type="text" value="50.00"/>

Balance Change	
Test Load	Indication
0	<input type="text" value="0.00"/>

Remarks

I copied the comments in red displayed at the top to remarks and hit Save.

All subsequent remarks are saved the same way.

It show who and when each remarks entry is made.

Remarks

7/2/2008 1:25:04 PM : Mary Vick : Indication for increasing test load 100 outside of tolerance range of +/- 0.020 Checktest submitted and DISAPPROVED

Previous

Save

Cancel

Void

Submit

Equipment Status Report for Lab Scales

Indicates Next Test Date

Equipment Status Reports Parameters for Report Generation

Equipment Type	Serial #	Make	Model	TRN	Test Date	Status	Result	Next Test Date
244120	Kansas City				MO			
Laboratory Scale	1175	Ohaus	GT4400-G	14028	6/27/2008	Closed	Pass	12/24/2008
244121	KANSAS CITY				MO			
Laboratory Scale	602	Ohaus	GT4100G	14027	6/27/2008	Closed	Fail	12/24/2008
Laboratory Scale	G-31359	Mettler	PM400/49	14043	7/1/2008	Closed	Pass	12/28/2008
400010	Aberdeen				SD			
Laboratory Scale	1125322606P	Ohaus	ECD112	14022	6/20/2008	Closed	Fail	12/17/2008
Laboratory Scale	1009	Ohaus	GT400E	14021	6/18/2008	Closed	Pass	12/15/2008

Redo Generate

More Examples

➤ *PM 2200*

➤ *PJ400*

FGIS-904 Laboratory Scale Test Test Information

Test Reference Number: 14053

Checktest submitted and APPROVED

1. **Class F weights** used in this exam were last tested on this date

2. **Zero balance** the scale before the sensitivity test, shift test and increasing load test.
 Yes No

3. **Sensitivity** Mechanical scales only. At zero, add weight equal to one scale division. The scale indicator should move to the top of the trig loop (or readings should change by at least one division). Repeat the test at the scale's maximum test load.
 Yes No

Shift Test				
Shift Load Test: <input type="text" value="1050.0"/> expected shift test reading)				
	Quadrant 1		Quadrant 2	
Reading	<input type="text" value="1050.0"/>	<input type="text" value="1050.0"/>	<input type="text" value="1050.0"/>	<input type="text" value="1050.0"/>

Increasing Load							
Test Load	Indication	Test Load	Indication	Test Load	Indication	Test Load	Indication
1g	<input type="text" value="1.00"/>	7g	<input type="text" value="7.00"/>	100g	<input type="text" value="100.00"/>	500	<input type="text" value="500.0"/>
2	<input type="text" value="2.00"/>	8	<input type="text" value="8.00"/>	150	<input type="text" value="150.00"/>	600	<input type="text" value="600.0"/>
3	<input type="text" value="3.00"/>	9	<input type="text" value="9.00"/>	200	<input type="text" value="200.0"/>	700	<input type="text" value="700.1"/>
4	<input type="text" value="4.00"/>	10	<input type="text" value="10.00"/>	250	<input type="text" value="250.0"/>	1000	<input type="text" value="1000.0"/>
5	<input type="text" value="5.00"/>	20	<input type="text" value="20.00"/>	300	<input type="text" value="300.0"/>	1500	<input type="text" value="1500.1"/>
6	<input type="text" value="6.00"/>	50	<input type="text" value="50.00"/>	400	<input type="text" value="400.0"/>	2100	<input type="text" value="2099.9"/>

Decreasing Load	
Test Load	Indication
400	<input type="text" value="400.0"/>
200	<input type="text" value="200.0"/>
50	<input type="text" value="50.00"/>

Balance Change	
Test Load	Indication
0	<input type="text" value="0.00"/>

Approved Disapproved
Tested By:

Remarks

PM2200

Dual Range Class II

FGIS-904 Laboratory Scale Test General Information

Test Reference Number: 14053

Test Date:

Service Point:

FO: FGIS - Wichita Field Office OA LOCATION: KANSAS CITY, MO

Laboratory Scale

SP Code	Serial #	Make	Model
141341	K98331	Mettler	PM2200-1/49

Capacity: g X Division Size 1 (d): g

Dual Range Scales Only:

Break Point: g Division Size 2 (d): g

Class II

Class II		
No. of Div.	Tol.	
0-5,000	1d	
5,001-20,000	2d	
20,001+	3d	

 Class III

Class III		
No. of Div.	Tol.	
0-500	1d	
501-2,000	2d	
2,001-4,000	3d	
4,001+	5d	

Type	d	Test Loads Required	Tol.
<input checked="" type="radio"/> Precision, etc.	0.01 g	1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 20, 50, 100 g	0.02 g
<input type="radio"/> Moisture, etc.	0.1	above plus 150, 200, 250, 300, 400 g	0.2g
<input type="radio"/> General, etc.	1	above plus 500, 600, 700, 1000, 1500 g	1



PJ400

**FGIS-904 Laboratory Scale Test
Test Information**

General Test

Test Reference Number: 14055

Checktest submitted and APPROVED

1. **Class F weights** used in this exam were last tested on this date

2. **Zero balance** the scale before the sensitivity test, shift test and increasing load test.
 Yes No

3. **Sensitivity** Mechanical scales only. At zero, add weight equal to one scale division. The scale indicator should move to the top of the trig loop (or readings should change by at least one division). Repeat the test at the scale's maximum test load.
 Yes No N/A

You are here: ECT Menu > FGIS-904 Laboratory Scale Form

**FGIS-904 Laboratory Scale Test
General Information**

General Test

Test Reference Number: 14055

Test Date:

Service Point

FO OA LOCATION
 FGIS - Wichita Field Office KANSAS CITY, MO

Laboratory Scale

SP Code	Serial #	Make	Model
461590	H-43095	Mettler	PJ400/49

Capacity g X Division Size 1 (d) g

Dual Range Scales Only:
 Break Point g Division Size 2 (d) g

<input checked="" type="radio"/> Class II	<table border="1"> <thead> <tr> <th colspan="2">Class II</th> </tr> <tr> <th>No. of Div.</th> <th>Tol.</th> </tr> </thead> <tbody> <tr> <td>0-5,000</td> <td>1d</td> </tr> <tr> <td>5,001-20,000</td> <td>2d</td> </tr> <tr> <td>20,001+</td> <td>3d</td> </tr> </tbody> </table>	Class II		No. of Div.	Tol.	0-5,000	1d	5,001-20,000	2d	20,001+	3d	<input type="radio"/> Class III	<table border="1"> <thead> <tr> <th colspan="2">Class III</th> </tr> <tr> <th>No. of Div.</th> <th>Tol.</th> </tr> </thead> <tbody> <tr> <td>0-500</td> <td>1d</td> </tr> <tr> <td>501-2,000</td> <td>2d</td> </tr> <tr> <td>2,001-4,000</td> <td>3d</td> </tr> <tr> <td>4,001+</td> <td>5d</td> </tr> </tbody> </table>	Class III		No. of Div.	Tol.	0-500	1d	501-2,000	2d	2,001-4,000	3d	4,001+	5d
Class II																									
No. of Div.	Tol.																								
0-5,000	1d																								
5,001-20,000	2d																								
20,001+	3d																								
Class III																									
No. of Div.	Tol.																								
0-500	1d																								
501-2,000	2d																								
2,001-4,000	3d																								
4,001+	5d																								

Type	d	Test Loads Required	Tol.
<input type="radio"/> Precision, etc.	0.01 g	1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 20, 50, 100 g	0.02 g
<input type="radio"/> Moisture, etc.	0.1	above plus 150, 200, 250, 300, 400 g	0.2g
<input type="radio"/> General, etc.	1	above plus 500, 600, 700, 1000, 1500 g	1

Save Cancel Void Next

Shift Test				
Shift Load Test:	<input type="text" value="200.00"/>	(used as expected shift test reading)		
	Quadrant 1	Quadrant 2	Quadrant 3	Quadrant 4
Reading	<input type="text" value="199.98"/>	<input type="text" value="199.98"/>	<input type="text" value="199.98"/>	<input type="text" value="199.98"/>

Increasing Load							
Test Load	Indication	Test Load	Indication	Test Load	Indication	Test Load	Indication
1g	<input type="text" value="1.00"/>	7g	<input type="text" value="7.00"/>	100g	<input type="text" value="99.99"/>	500	<input type="text" value=""/>
2	<input type="text" value="2.00"/>	8	<input type="text" value="8.00"/>	150	<input type="text" value="149.99"/>	600	<input type="text" value=""/>
3	<input type="text" value="3.00"/>	9	<input type="text" value="9.00"/>	200	<input type="text" value="199.98"/>	700	<input type="text" value=""/>
4	<input type="text" value="4.00"/>	10	<input type="text" value="10.00"/>	250	<input type="text" value="249.98"/>	1000	<input type="text" value=""/>
5	<input type="text" value="5.00"/>	20	<input type="text" value="20.00"/>	300	<input type="text" value="299.97"/>	1500	<input type="text" value=""/>
6	<input type="text" value="6.00"/>	50	<input type="text" value="50.00"/>	400	<input type="text" value="399.97"/>		<input type="text" value=""/>

Decreasing Load	
Test Load	Indication
400	<input type="text" value="399.97"/>
200	<input type="text" value="199.98"/>
50	<input type="text" value="50.00"/>

Balance Change	
Test Load	Indication
0	<input type="text" value="0.00"/>

Approved Disapproved Tested By: